ABOUT US


Xiamen Wiz Biotech Co., Ltd. is a high-tech biomedical company, devoted to the field of rapid diagnostic reagents and instrument. Wiz is located in Xiamen, China, established in July 2013, which is a listed enterprises in National Equities and Quotations. As a technology-driven company that prides itself on its scientific excellence, Wiz focused on technological innovation and product innovation with 15 authorized patents, each of our products embodies the creativity of our excellent scientists, who are working hard to continuously bring novel products to the China and world markets. SARS-CoV-2 rapid detection series of products, including SARS-CoV-2 antigen, antibody and influenza differential detection, are produced by Wizbiotech, which can easily and quickly carry out screening of SARS-CoV-2 in a large number of people, adding strength to epidemic prevention and control.

Xiamen Wiz Biotech Co., Ltd.

COVID-19 Antigen Self-Testing Kit

Self-testing


SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

- ✓ Anterior nasal swab specimen, noninvasive
- ✓ Very simple to use
- ✓ Convenient, no instrument required
- ✓ Rapid, results within 15~20 minutes
- ✓ Cost-efficient

<http://en.wizbiotech.com/>

CE 1434

QUICK REFERENCE GUIDE

CLINICAL PERFORMANCE

Test Results	Reference PCR Results		
	Positive	Negative	Total
Positive	113	0	113
Negative	2	456	458
Total	115	456	571

Sensitivity: **98.26%** (95% C.I. 93.86%~99.79%)
 Positive Predictive Value: **100%** (95% C.I. 96.79%~100.00%)
 Overall Percent Agreement: **99.65%** (95% C.I. 98.74~99.96%)

Specificity: **100.00%** (95% C.I. 99.19%~100.00%)
 Negativity Predictive Value: **99.56%** (95% C.I. 98.43%~99.95%)

PACKING INFORMATION

Specifcaton	Box size(mm)	Carton size(mm)	Carton G.W.(kg)	Boxes/Carton	Tests/Carton
1Test	150*70*15mm	615*315*370mm	14kg	400	400
2Tests	150*65*30mm	675*320*320mm	11kg	200	400
5Tests	170*80*40mm	420*360*340mm	9kg	80	400
10Tests	150*100*80mm	420*320*420mm	7kg	40	400
20Tests	210*140*70mm	440*300*370mm	7kg	20	400

CERTIFICATE

CERTIFICATE

EC Certificate No. 1434-IVDD-492/2021

EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies
that manufactured by:

Xiamen Wiz Biotech Co., Ltd.,
3-4 Floor, NO.16 Building, Bio-medical Workshop,
2030 Wengjiao Xi Road, Haicang District, Xiamen City,
Fujian Province, 361026, P.R. China

in vitro diagnostic medical devices
for self-testing

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

51332801, 51332802, 51332803, 51332804

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 22.11.2021 to 27.05.2024

The date of issue of the Certificate: 22.11.2021

The date of the first issue of the Certificate: 22.11.2021

Issued under the Contract No. MD-77/2021
Application No: 130/2021
Certificate bears the qualified signature.
Warsaw, 22/11/2021
Module A1

Vice-President

POLISH CENTRE FOR TESTING AND CERTIFICATION 02-844 Warsaw, 469 Pulawska Street, tel. +48 22 46 45 200, e-mail: pcbc@pcbc.gov.pl

Certificate CN1942110

The management system of

Xiamen Wiz Biotech Co., Ltd.

3-4 Floor, NO. 16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road,
Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

has been assessed and certified as meeting the requirements of

ISO 13485:2016
EN ISO 13485:2016

For the following activities

The scope of registration appears on page 2 of this certificate.

This certificate is valid from 26 August 2020 until 3 September 2022
and remains valid subject to satisfactory surveillance audits.

Re certification audit due before 26 July 2022
Issue 2. Certified since 4 September 2019

This is a multi-site certification.
Additional site details are listed on the subsequent page.

Authorized by

SGS United Kingdom Ltd
Rosemore Business Park, Ebbw Vale, Gwent, NP23 5EH, UK
T +44 (0)151 350-6666 F +44 (0)151 350-6669 www.sgs.com

HC SGS 13485 2016 0118 M2

Page 1 of 2

This document is issued by the Company subject to the General Conditions of Certification. The user is responsible for the use of the certificate and for the validity of the information provided. The user is also responsible for the use of the certificate and for the validity of the information provided. The user is also responsible for the use of the certificate and for the validity of the information provided.

Certificate CN1942110, continued

Xiamen Wiz Biotech Co., Ltd.

ISO 13485:2016
EN ISO 13485:2016

Issue 2

Detailed scope

Design and manufacture of *in vitro* diagnostic medical devices
including hormone marker test kit, kidney disease marker test kit,
tumor marker test kit, infectious disease marker test kit,
gastrointestinal inflammation marker test kit, cardiac marker test kit,
immunochromatography analyzer.

Additional facilities

4 Floor, NO. B12 Building, Bio-medical Workshop, 2072 Wengjiao Xi Road,
Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

This document is issued by the Company subject to the General Conditions of Certification. The user is responsible for the use of the certificate and for the validity of the information provided. The user is also responsible for the use of the certificate and for the validity of the information provided. The user is also responsible for the use of the certificate and for the validity of the information provided.

Page 2 of 2

FACTORY

Our one hundred thousand grade clean and dust-free workshop is operated strictly under ISO 13485:2016 and GMP guidelines. With the automatic production line, efficient production process, strict quality control, we always produce the products timely to meet customers' request.

INSTRUCTIONS FOR USE

EN

EN

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

Instructions for Use

—For use at home self-test or Non-professional
—For use with nasal cavity (anterior nasal) swab specimen
—For In Vitro Diagnostic Use Only

INTENDED USE

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) is intended for the qualitative detection of SARS-CoV-2 Antigen (Nucleocapsid protein) which is in nasal cavity (anterior nasal) swab specimen from individuals with suspect COVID-19 infection. The test kit is intended for self-test or home test.

SUMMARY

COVID-19 is an acute respiratory infectious disease; People are prone to infection generally. Currently, the main cause of COVID-19 infection is contact with someone who is already infected with the SARS-CoV-2, and asymptomatic infected people may also transmit the virus. Studies have shown that symptoms of infection generally appear within 14 days, with most occurring within 3 to 7 days after infection. The main symptoms are fever, fatigue, loss of smell or taste, and a dry cough. In some cases, a stuffy nose, runny nose, muscle pain and diarrhea can also occur.

PRINCIPLE THE DETECTION

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) adopts immune lateral chromatography technology. When the test specimen contains SARS-CoV-2, the SARS-CoV-2 Antigen will react with the antibody coated on the test line (T) to appear a red band in the test line (T) area; when the content of the SARS-CoV-2 in the test specimen is too low or does not exist, the test line (T) area does not appear a red band. Regardless of whether the test specimen contains the SARS-CoV-2, a red band will appear in the area of the quality control line (C), which is the basis for judging whether the test is effective.

MAIN KIT COMPONENTS

Catalogue number	51332801	51332802	51332803	51332804
Specification	1 Test/kit	2 Tests/kit	5 Tests/kit	10 Tests/kit
Test device	1	2	5	10
Extraction tube	1	2	5	10
Disposable swab	1	2	5	10
Biohazard Waste Bag	1	2	5	10
Instructions for Use	1	1	1	1

Materials required but not provided: Timer

PRECAUTIONS

- For In Vitro Diagnostic Use Only.
- For use with nasal cavity (anterior nasal) swab specimen.
- For the detection of protein from SARS-CoV-2 only, not for detect other viruses or pathogens.
- The person under 18 years of age should be tested with the assistance of the legal guardian or authorized person.
- Keep the test kit or kit components out of the reach of children and pets before and after use.
- Test card packaging contains desiccant, it is forbidden to eat.
- The sample extraction solution in extraction tube contains chemical components. Direct contact should be avoided and eat prohibited. If the solution contacts the skin or mucosa or eye, please flush with copious amounts of water. Please contact your family doctor or professional or seek medical advice if necessary.
- Use of gloves and other protective equipment are recommended when conducting testing.
- The test kits should be stored according to storage conditions required in instruction for use. It is forbidden to use the test kit that is not stored as required.
- Don't use the test kit beyond its expiration date.
- Don't use kit components that have been opened or changed.
- Leave test card sealed in its foil pouch until just before use. Don't use if pouch is damaged or open.
- Disposable swabs are sterile products. Don't use if the swab packaging is damaged or opened.
- The use of the swab should strictly follow instruction for use; otherwise it may cause nasal cavity bleeding, swab rupture and retention and other risks.
- Don't immerse the swab in the extraction solution or other liquid before collect sample by Disposable swabs.
- Don't touch swab soft tip when handling the swab sample.
- Proper sample collection and handling are crucial for correct results.
- Don't mix components from different kit lots.

19. All kit components are single use items. Don't use with multiple specimens. Don't reuse the used test kit or kit components.

20. Before deciding to implement relevant treatment or management decisions, it is recommended to conduct communication with family doctors or professionals. Don't take medicine privately or any action that will endanger yourself or others.

21. The used test kit components and samples can be put into plastic bags together with ordinary household waste. If the test result is positive, you should carefully dispose the relevant waste components and samples, and thoroughly clean and disinfect the working surface to ensure hygiene. If there are special requires on waste in local laws and regulations, you must strictly observe it.

22. In view of the global epidemic of COVID-19, all actions should comply with the current measures and regulations of your country/region, to implement prevention and control measures scientifically, and protect yourself and others effectively.

23. Do not eat, drink or smoke in the area where handling specimens or test kits.

STORAGE CONDITIONS AND SHELF LIFE

The test kit should be stored conditions of 2°C~30°C, dry and out of direct sunlight (Don't freeze the kit or its components).

The shelf life of the kit is 12 months.

The test card should be used within 60 minutes after opening the aluminum foil bag.

For the kit expiration date, please refer to the product label.

FREQUENTLY ASKED QUESTIONS (FAQ)

1. What are the known or potential benefits of product testing?

— The test results can help your family doctor or professional make accurate or effective recommendations.

— Test results may help limit the spread of COVID-19 to your family and others in your community.

2. What are the known or potential risks in product testing?

— Discomforts may occur during specimen collection.

— Incorrect test results may be obtained.

3. When should/can I test myself?

You can test yourself when you have suspected symptoms of COVID-19. Studies have shown that COVID-19 infected person have a high virus load in the first four days of illness, making it easier to detect.

4. What's the difference between Antigen and Molecular test?

Currently, there are several SARS-CoV-2 test methods. Molecular test (also known as PCR test) detect the genetic material of the virus, and Antigen test detect the proteins in the virus.

5. What factors will affect the test results? What should I pay attention to?

— For anterior nasal swab specimen only.

— The sample must not contain bubbles when dripping.

— Do not add too much or too little sample.

— Test immediately after sample collection.

— Strictly follow the instructions for use.

6. No red line band on test card or abnormal fluid flow? What is the reason?

It should be clear that the test result is invalid. The reasons are as follows:

— The table on which the test card is placed is not smooth, affecting the flow of liquid.

— Drop sample volume does not meet the requirements specified in the instructions for use.

— The test card is damp.

7. I have taken the test, but I don't see the control line (C). What should I do?

Your test result is invalid. Please repeat the test strictly according to the instructions for use.

8. Unsure about the test result, what should I do?

For uncertain results can be retested. If you're still unsure of the test result, please contact the nearest medical institution according to the advice of your local government.

9. If result is positive, what should I do?

If your test result is positive, you may be infected with COVID-19. You should take the necessary measures (such as quarantine, report, retest, etc.) as required by the local government, and contact the nearest medical institution for next action.

10. If result is negative, what should I do?

If the test kit only clearly shows the control line band, this may mean that the test result is negative or that the virus is too low to be detected. If you still have COVID-19 symptoms (headache, fever, loss of smell or taste, etc.), please consult your family doctor or medical institution recommended by local government.

POSITIVE VALUE

Positive value/limit of detection: 1.7×10^4 TCID₅₀/mL

Select the confirmed inactivated SARS-CoV-2 medium (concentration 3.4×10^4 TCID₅₀/mL), and use gradient dilution method to find out the virus medium which is reach the critical value of the detection. That is repeating the test for 20 times and the test result is positive for at least 19 times.

CLINICAL PERFORMANCE

The clinical performance characteristic of the SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) was established with 571 anterior nasal swab specimens. Specimens were collected and tested according to the requirements of the Instructions for Use. The storage, transportation and detection of samples after collection met the relevant requirements of the Instructions for Use. At the same time, the nasopharyngeal swab specimens of the same donor were detected by nucleic acid detection reagent (Vitalassay qPCR SARS-CoV-2). WIZ'S SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) was compared with reference reagent PCR results.

Clinical performance of the WIZ'S SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) compared to RT-PCR

WIZ Results	Reference PCR Results		
	Positive	Negative	Total
Positive	113	0	113
Negative	2	456	458
Total	115	456	571

Sensitivity: 98.26% (95%CI: 93.86%~99.79%)

Specificity: 100.00% (95%CI: 99.19%~100.00%)

Positive Predictive Value: 100% (95%CI: 96.79%~100.00%)

Negativity Predictive Value: 99.56% (95%CI: 98.43%~99.95%)

Overall Percent Agreement: 99.65% (95%CI: 98.74~99.96%)

PERFORMANCE CHARACTERISTICS

- Using enterprise reference for testing, the results meet the requirements of enterprise reference.
- Cross reaction

Name (Microorganism)	Concentration	Test result
Influenza B/Yamagata	1.83×10^4 TCID ₅₀ /mL	Negative
Influenza B/Victoria	2.07×10^4 TCID ₅₀ /mL	Negative
Influenza A H3N1	1.00×10^4 TCID ₅₀ /mL	Negative
Influenza A H3N2	1.15×10^4 TCID ₅₀ /mL	Negative
Influenza A H5N1	1.32×10^4 TCID ₅₀ /mL	Negative
Avian Influenza A (H7N9)	1.60×10^4 TCID ₅₀ /mL	Negative
SARS Coronavirus	2.14×10^4 TCID ₅₀ /mL	Negative
Adenovirus 1	1.39×10^4 TCID ₅₀ /mL	Negative
Adenovirus 3	1.24×10^4 TCID ₅₀ /mL	Negative
Adenovirus 7	1.87×10^4 TCID ₅₀ /mL	Negative
Human coronavirus 229E	2.00×10^4 TCID ₅₀ /mL	Negative
Human coronavirus OC43	2.34×10^4 TCID ₅₀ /mL	Negative
Human coronavirus NL63	2.00×10^4 TCID ₅₀ /mL	Negative
Human coronavirus HKU1	2.00×10^4 TCID ₅₀ /mL	Negative
MERS-coronavirus	1.00×10^4 TCID ₅₀ /mL	Negative
Cytomegalovirus	1.00×10^4 TCID ₅₀ /mL	Negative
Enterovirus 71	2.55×10^4 TCID ₅₀ /mL	Negative
Human Parainfluenza virus 1	1.35×10^4 TCID ₅₀ /mL	Negative
Human Parainfluenza virus 2	6.31×10^4 TCID ₅₀ /mL	Negative
Human Parainfluenza virus 3	3.25×10^4 TCID ₅₀ /mL	Negative
Human Parainfluenza virus 4	3.31×10^4 TCID ₅₀ /mL	Negative
Measles virus	6.31×10^4 TCID ₅₀ /mL	Negative
Mumps virus	6.31×10^4 TCID ₅₀ /mL	Negative
Respiratory syncytial virus	2.00×10^4 TCID ₅₀ /mL	Negative
Rhinovirus 1A	1.26×10^4 TCID ₅₀ /mL	Negative
Rhinovirus	1.20×10^4 TCID ₅₀ /mL	Negative
Epstein Barr Virus	2.18×10^4 TCID ₅₀ /mL	Negative
Varicella zoster virus	1.00×10^4 TCID ₅₀ /mL	Negative

Bacillus pertussis	1.30×10^4 CFU/mL	Negative
Chlamydia pneumoniae	1.00×10^4 CFU/mL	Negative
Escherichia coli	1.00×10^4 CFU/mL	Negative
Haemophilus influenzae	1.20×10^4 CFU/mL	Negative
Mycobacterium binding	1.00×10^4 CFU/mL	Negative
Mycoplasma Pneumoniae	1.00×10^4 CFU/mL	Negative
Candida Albicans	1.00×10^4 CFU/mL	Negative
Neisseria meningitidis	1.00×10^4 CFU/mL	Negative
Neisseria gonorrhoeae	1.00×10^4 CFU/mL	Negative
Pseudomonas aeruginosa	1.70×10^4 CFU/mL	Negative
Staphylococcus aureus	2.20×10^4 CFU/mL	Negative
Streptococcus pneumoniae	1.00×10^4 CFU/mL	Negative
Streptococcus pyogenes	1.28×10^4 CFU/mL	Negative
Streptococcus salivarius	1.00×10^4 CFU/mL	Negative
Legionella Pneumophila	1.98×10^4 CFU/mL	Negative

Interference Substances

Interference Substances	Concentration	Negative interference results	Positive interference results
Mucin	5%	Negative	Positive
Whole blood	5% (V/V)	Negative	Positive
α -interferon	500 thousand IU/mL	Negative	Positive
Zanamivir	500mg/mL	Negative	Positive
Ribavirin	20ug/mL	Negative	Positive
Osetamivir	5ug/mL	Negative	Positive
Peramivir	0.2mg/mL	Negative	Positive
Lopinavir	8 mg/mL	Negative	Positive
Ritonavir	550ug/mL	Negative	Positive
Umifenovir	4ug/mL	Negative	Positive
Levofloxacin	30ug/mL	Negative	Positive
Azithromycin	4.5ug/mL	Negative	Positive
Ceftriaxone	0.8 mg/mL	Negative	Positive
Meropenem	1.3mg/mL	Negative	Positive
Tobramycin	4mg/mL	Negative	Positive
Phenylephrine	20ug/mL	Negative	Positive
Glymetazoline	0.1mg/mL	Negative	Positive
Beclomethasone	0.1mg/mL	Negative	Positive
Dexamethasone	2 mg/mL	Negative	Positive
Fluticasone	0.1mg/mL	Negative	Positive
Triamcinolone acetonide	10.5mg/mL	Negative	Positive
Budesonide	2.75mg/mL	Negative	Positive
Mometasone	10mg/mL	Negative	Positive
Fluticasone	55ug/mL	Negative	Positive
Histamine Hydrochloride	10mg/mL	Negative	Positive
Sodium chloride	5%	Negative	Positive

Hook effect

Within the concentration of 3.4×10^4 TCID₅₀/mL for cell culture medium of SARS-CoV-2 Antigen, the test results of this product showed no Hook effect.

Repeatability

The product was established by using internal reference, there were no differences observed within-run, between-run.

LITERATURE REFERENCES

- NMPA. The Technical Key Points for Coronavirus (COVID-19) Antigen-antibody Detection Reagent Registration Review (Trial). (2020).
- Xu Zhao, Li Ran. Analysis on the Risk Management of in Vitro Diagnostic Reagents[J]. China Medical Device Information. 2020, 26(13):8-10.
- Wu Jinhui, Meng Li. Immunocolloidal Gold Technology: Advances and Application[J]. Chinese Agricultural Science Bulletin. 2019, 35(13): 146-151.
- Li Yorguin, Yang Ruffu. Rapid test of immunocolloidal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2003, 31(1): 74-78.

TEST PROCEDURE

It is important to read the Instructions for Use carefully and follow the steps in the correct order. It takes about 15 minutes to prepare before each test, and the results can be obtained after 30 minutes.

Test Procedure	Preparatory work	Sample collecting	Sample processing	Sample testing
Chapter	I	II	III	IV

I Preparatory work

- Please use the test kit at room temperature (15°C~30°C). If the test kit was previously stored in a cool place (temperature below 15°C), please balance at 15~30°C for 30 minutes before test.
- Prepare a Timer (such as watch, clock), tissue, hand sanitizer/soap and warm water.
- Please read the Instruction for Use carefully.

5. Check the components in the box and verify that there is not any damage or breakage.

Test device (sealed package)	Disposable swab (sealed package)	Extraction tube (filled with extraction solution)	Biohazard Waste Bag	Timer (prepare by yourself)

Warning!

If the sealed package is damaged, don't use it. Please replace it with a new kit. If you have cough symptoms, get tested in private.


II Sample collecting

- Take out the sample extraction tube; unscrew the tube cap of the extraction tube.


- Place the extraction tube on the bracket (attached to the box) to avoid liquid spillage.

- Tear off the swab's packing from stick end and take out the swab.


- As shown in the picture, wipe two nostrils with a swab:
(1) Insert swab soft head into nostril less than 1 inch (Normally about 0.5~0.75 inch).
(2) Gently rotate the swab against the nasal wall with moderate force at least five times.
(3) Use the same swab to repeat the collection procedure in another nostril.


III Sample processing

- Take out the extraction tube, dip the swab soft head into the extraction solution and immerse it in the liquid.

- Press the swab soft head tightly against the inner wall of the extraction tube, and rotate the swab clockwise or anticlockwise for about 10 times.


- Squeeze the swab head along the inner wall of the sample extraction tube under the liquid in the tube as much as possible, take out the swab.

IV Sample testing

- Tear off the aluminum foil bag, take out the test card and place it horizontally on the test desk.


Warning!

The platform should be in a horizontal and stable state, and tilt and shake are strictly prohibited.

- Unplug the adding sample hole cover of the extraction tube.


- Gently squeeze the extraction tube, and drop 2 drops liquid vertically into the sample well of the test card.

Warning!

The existence of bubbles in the extraction tube may lead to the wrong sample volume and inaccurate test results. If there are bubbles in the extraction tube, gently shake the extraction tube to squeeze out part of the liquid so as to remove the bubbles.

- Start timing, read the test results at 15 minutes. Do not read result before 15 minutes or after 30 minutes.

- After test is completed, put all test kit materials into the biohazard waste bag and dispose it according to the local biohazard waste disposal policy.

- Rewash hands thoroughly (at least 20 seconds) with soap and warm water/hand sanitizer.

INTERPRETATION OF THE RESULTS

Negative: The quality control line (C line) appears red band, while the test line (T line) doesn't appear red band.	
Positive: Both the quality control line (C line) and the test line (T line) appear red bands.	
Caution! Regardless of the shades of color within the specified detection time, the result should be judged as positive result.	

Invalid: Once the red line in the control line (C line) region doesn't appear which will be treated as invalid.	
Warning! Invalid results may indicate that you did not strictly follow the operation steps of instruction for use. Please read the instruction for use carefully again and select a new kit for testing again.	

PROCESSING OF TEST RESULTS

- If the test result is positive:
 - COVID-19 infection is present currently.
 - Contact your doctor or local health department immediately.
 - Comply with the self-quarantine requirement and protection guidelines in your area.
 - Carry out PCR test for confirmation.
- If the test result is negative:
 - You need to continue to take measures related to contact with others and self-protection.
 - Even if your test result is negative, it is also possible that there is an infection.
 - If you still suspect, please repeat the test after 1-2 days. On account of the coronavirus cannot be accurately detected at every stages of infection.
- If the test result is invalid:
 - It may be caused by incorrect operation in the detection process.
 - Please repeat the test.
 - If the test result is still invalid, please contact your doctor or COVID-19 testing center.

QUALITY CONTROL

The quality control line is a key point of test kit and is used to control the procedure. The quality control line appears, which indicates that the test has been performed correctly and the test kit has reacted.

SYMBOLS

Symbol	Used for	Symbol	Used for	Symbol	Used for
	Consult instructions for use		Test kit		Manufacturer
	In Vitro Diagnostic Medical Device		Use by date		Do not use
	Store at 2°C – 8°C		Lot/unique number		Batch code
	Authorized Representative of the European Commission		Location		Keep dry
	Use the product about the package is damaged		Biological risks		Test devices
	In Vitro Diagnostics		Disposable results		Individual Waste

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RIFERIMENTI LETTERARI

- Wu, The Technical Key Points for Coronavirus (COVID-19) Antigen-Antibody Detection Reagent Registration Review (Trial), (2020).
- Xu Chao, U Ren. Analysis on the Risk Management of in Vitro Diagnostic Reagents[J]. China Medical Device Information, 2020, 26(13): 8-10.
- Wu Anhu, Meng Li. Immunocolloidal Gold Technology: Advances and Application[J]. Chinese Agricultural Science Bulletin, 2019, 35(13): 146-151.
- Li Yongping, Yang Ruliu. Rapid test of immunocolloidal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2003, 23(1): 74-78.

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REFERENCIAS DE LITERATURA

- NMFA. The Technical Key Points for Coronavirus (COVID-19) Antigen-antibody Detection Reagent Registration Review (Trial) [2020].
- Xu Chao, Li Ran. Analysis on the Risk Management of in Vitro Diagnostic Reagents[J]. China Medical Device Information, 2020, 26(12):75-79.
- Wu Jinhui, Meng Li. Immunocolloidal Gold Technology: Advances and Application[J]. Chinese Agricultural Science Bulletin, 2019, 35(13): 146-151.
- Li Yonggen, Tang Ruli. Rapid test of immunocolloidal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2003, 31(1): 74-78.

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LITERATUR REFERENZEN

- NMPA. The Technical Key Points for Coronavirus (COVID-19) Antigen-antibody Detection Reagent Registration Review (Trial) (2020).
- Xu Chao, Li Ren. Analysis on the Risk Management of In Vitro Diagnostic Reagents[J]. China Medical Device Information, 2020, 24(13):8-10.
- Wu Jinhua, Meng Li. Immunocolloidal Gold Technology: Advances and Application[J]. Chinese Agricultural Science Bulletin, 2019, 35(13): 146-151.
- Li Yongqiang, Yang Rukun. Rapid test of immunocolloidal gold with membrane as solid phase carrier[J]. Progress in Microbiology and Immunology, 2000, 33(1): 74-76.

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Xiamen WIZ Biotech Co., Ltd.
Address: 3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China
Tel: +86-592-6808278 Fax: +86-592-6808279

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0123 (According to Directive 93/42/EEC)

Address: No. 10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

EC REP Name: Shanghai International Holding Corp. GmbH(Europe)
EC REP Address: Eiffelstrasse 80, 20537 Hamburg, Germany

CE-01.0301

CONTACT US

CONTACT US

XIAMEN WIZ BIOTECH CO.,LTD

Tel:

+86-592-6808278

Fax:

+86-592-6808279

E-mail:

sales@wizbiotech.com

Address:

3-4 Floor, NO.16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

www.wizbiotech.com

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Zhejiang Gongdong Medical Technology Co., Ltd.