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Data: 04/12/2021

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Ct-value	No. of Samples	No. of true positive Rapid Test Samples	Sensitivity of SARS-CoV-2 Antigen Rapid Test
≤30	82	82	100%
≤32	94	92	97.9%
≤34	102	98	96.1%
≤36	109	103	94.5%

2. Analytical Results with correlation to Ct-values of the negative samples:

No. of Samples	No. of true negative Rapid Test Samples	Sensitivity of SARS-CoV-2 Antigen Rapid Test
300	300	100%

- France:**
SPIRAL Evaluation with good results: Sensitivity 97.1%, Specificity 100%
- Malaysia:**
IMR(Institute for Medical Research) Evaluation with good results: Sensitivity 96.0%, Specificity 100%
- Thailand:**
Evaluation with good results: Sensitivity 100%, Specificity 100%

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No. Q5 095295 0001 Rev. 00

Holder of Certificate:

Hangzhou Sejoy Electronics & Instruments Co., Ltd.
Area C, Building 2, No. 365, Wuzhou Road
Yuhang Economic Development Zone
311100 Hangzhou City, Zhejiang
PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate:

Design and Development, Production and Distribution of In Vitro Diagnostic Medical Device based on Immunochromatography, Dry Chemistry and Electrochemistry Method, Include Instrument, Test Strip and Control Solution

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:Q5 095295 0001 Rev. 00

Report No.:

SH20167601

Valid from:

2020-10-30

Valid until:

2023-10-29

Date,

2020-10-30

Christoph Dicks

Head of Certification/Notified Body

CERTIFICATE

EC Certificate No. 1434-IVDD-474/2021

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

**Hangzhou Sejoy Electronics & Instruments Co., Ltd
Area C, Building 2, No. 365, Wuzhou Road, Yuhang Economic
Development Zone, 311100 Hangzhou City, Zhejiang, China**

in vitro diagnostic medical devices
for self-testing

**SARS-CoV-2 Antigen Rapid Test Cassette
COVG-602ST**

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 22.10.2021 to 27.05.2024

The date of issue of the Certificate: 22.10.2021

The date of the first issue of the Certificate: 22.10.2021

Issued under the Contract No. MD-100/2021
Application No: 192/2021
Certificate bears the qualified signature.
Warsaw, 22/10/2021
Module A1
FBM-30-E_10

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EU DECLARATION OF CONFORMITY

Manufacturer: Hangzhou Sejoy Electronics & Instruments Co.,Ltd.
Area C, Building 2, No.365, Wuzhou Road,
Yuhang Economic Development Zone,
311100 Hangzhou,Zhejiang,China

European Authorized Representative: Shanghai International Holding Corp.GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Product Name: SARS-CoV-2 Antigen Rapid Test Cassette

Model: COVG-602ST

Classification: Self-testing device not listed under Annex II of Directive 98/79/EC

Notified Body: Polish Centre for Testing and Certification
469 Puławska Street, 02-844 Warsaw

Notified Body No.: 1434

EC Certificate No.: 1434-IVDD-474/2021

Conformity assessment route: Annex III section 6 of Directive 98/79/EC

Applicable Standards: EN ISO 13485:2016, EN ISO 14971:2012,
EN 13532:2002, EN ISO 23640:2015, EN ISO 13612:2002,
EN ISO 17511:2003, EN 13975:2003,
EN ISO 18113-1:2011, EN ISO 18113-4:2011,
EN ISO 15223-1:2016, EN 13641:2002,EN 62366:2008

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

Hangzhou, October 23, 2021

Place, date

杭州世佳电子有限公司
HANGZHOU SEJOY ELECTRONICS & INSTRUMENTS CO.,LTD.

General Manager

Legally binding signature, Position

Statement on the monitoring of SARS-CoV-2 variants

The SARS-CoV-2 Antigen Rapid Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2 antigen in human nasal swab. The identification is based on the monoclonal antibodies specific for the Nucleocapsid (N) Protein of SARS-CoV-2. It is intended to aid in the rapid differential diagnosis of COVID-19 infection.

Recently, the SARS-CoV-2 has discovered the newest SARS-CoV-2 variant "Omicron", whose Pango lineage is B.1.1.529. The Sejoy urgently established a special verification team to monitor and analyze the genetic data of the newly discovered SARS-CoV-2 variant; The Peptide probe sequence comparison results of the marketed products confirmed that the SARS-CoV-2 Antigen Rapid Test Cassette (Ref.:COVG-602ST) that has been marketed by Hangzhou Sejoy Electronics & Instruments Co., Ltd. has no missed detection against the above-mentioned variant and still ensure the accuracy and sensitivity of the detection reagents.

Up to now, our company has monitored and analyzed the genetic data of major epidemic SARS-CoV-2 variants, including Alpha variant (B.1.1.7), Beta variant (B.1.351), Gamma variant (P.1) and Delta variant (B.1.617.2), Omicron variant (B.1.1.529), our company will continue to pay attention to the variant of the SARS-CoV-2 to ensure that our company's SARS-CoV-2 Antigen Rapid Test Cassette (Ref.:COVG-602ST) will not miss detection and ensure the sensitivity, accuracy and specificity are not affected.

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